# DSC Q Series 操作訓練

**2017.11.07**

**Allen Chang**

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DSC  
Q20/Q10

DSC  
Q200/Q100

DSC  
Q2000/Q1000

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QSeries - [Q2000-1584 - DSC Q2000@Mfg.dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 22.37°C

### Experiment

#### Standard Sequence

Sequence No. 12  
Run 1:

Summary Procedure Notes

Procedure Summary  
 Mode: Standard  
 Test: Custom

Sample Information  
 Sample Name: ITRI PET  
 Pan Type: Tzero Aluminum  
 Sample Size: 6.280 mg Pan No. 1 Ref. 1  
 Pan Mass: 0.000 mg (Sample) 0.000 mg (Reference)  
 Comments:

Data File Name: \\Ta-lab\ta\Data\DSC\LA WRENCE\TRIM\ITRI PET.001  
 Network Drive

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	22.37 °C
Heat Flow	0.173 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	21.25 °C
Heater Temperature	20.99 °C
Sample Purge Flow	50.01 mL/min

#	Running Segment Description
1	↕↑ Equilibrate at 0.00 °C
2	↕↘ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↕↘ Ramp 10.000 °C/min to 0.00 °C
5	↕↑ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↕↘ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 22.60 °C

### Experiment

Standard Sequence

Sequence No. 12  
Run 1:

Summary Procedure Notes

Procedure Information

Test: Custom

Notes:

Method Name: Editor... Advanced... Post Test...

#	Segment Description
1	↕ Equilibrate at 0.00 °C
2	↘ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↘ Ramp 10.000 °C/min to 0.00 °C
5	↕ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↘ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	22.60 °C
Heat Flow	0.164 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	21.69 °C
Heater Temperature	21.23 °C
Sample Purge Flow	49.99 mL/min

#	Running Segment Description
1	↕ Equilibrate at 0.00 °C
2	↘ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↘ Ramp 10.000 °C/min to 0.00 °C
5	↕ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↘ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3

### Advanced Parameters

Data Sampling Interval: 0.20 s/pt

Load Temperature Range 0.00 to 0.00 °C

OK Cancel Help

### Post Test Parameters

Method End Conditions

Unload Temperature Range 35.00 to 45.00 °C

Delay Time 0.00 min

Discard pan in waste bin

OK Cancel Help

### Method

Method Contents

Name: [ ]

#	Segment Description
1	↕ Equilibrate at 0.00 °C
2	↘ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↘ Ramp 10.000 °C/min to 0.00 °C
5	↕ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↘ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3

OK Cancel Help

Segment list

- ↕ Jump
- ↕ Equilibrate
- ↕ Initial temperature
- ↘ Ramp
- ↕ Isothermal
- ↕ Step
- ↕ Increment temperature
- ↕ Repeat
- ↕ Repeat until
- ↕ Abort next segment on limit
- ↕ Sampling interval

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Temperature (°C)

01 97.90 min. Append Apply Cancel Help

Ready

Completed Standard Seg 0 in Run 1 15:41:22

TA

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Start  
Stop  
Reject  
Hold  
Resume  
Lid  
Store  
Gas  
Event  
Air Cool  
PCA Shutter  
Go To Standby Temp  
Autosampler  
LNCS  
Shutdown Instrument

Complete Temp: 22.95°C

Summary Procedure Notes

Notes

Operator: RICO

Extended Text

Mass Flow Control Settings

Sample: #1 - Nitrogen Flow Rate: 50 mL/min

Auto Analysis

Autoanalyze

Analysis Macro

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	22.95 °C
Heat Flow	0.158 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	22.11 °C
Heater Temperature	21.59 °C
Sample Purge Flow	49.98 mL/min

#	Running Segment Description
1	↑ Equilibrate at 0.00 °C
2	↘ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↗ Ramp 10.000 °C/min to 0.00 °C
5	↑ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↘ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3

01 97.90 min Append Apply Cancel Help

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QSeries - [Q2000-1584 - DSC Q2000@Mfg.dsc]

Control Experimental Calibrate Tools View Window Help

Wizard... Ctrl + W  
Modify Running Method... Ctrl + M

Sequence File  
Sequence Run

Standard Sequence

Sequence No. 12  
Run 1

Notes

Operator: RICO

Extended Text

Mass Flow Control Settings

Sample: #1 - Nitrogen Flow Rate: 50 mL/min

Auto Analysis

Autoanalyze

Analysis Macro

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	23.21 °C
Heat Flow	0.155 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	22.36 °C
Heater Temperature	21.84 °C
Sample Purge Flow	49.98 mL/min

#	Running Segment Description
1	↑ Equilibrate at 0.00 °C
2	↓ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↑ Ramp 10.000 °C/min to 0.00 °C
5	↑ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↓ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3

01 97.90 min. Append Apply Cancel Help

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Calibration Wizard  
Analysis...  
AutoCalibration

Autosampler...  
Cell/Temperature Table...  
Indium Tracking  
Report...  
Mass Flow Control...  
Touchscreen

Experiment  
Standard Sequence

Sequence No. 12  
Run 1

Procedure Notes

RICO

Mass Flow Control Settings  
Sample #1 - Nitrogen Flow Rate 50 mL/min

Auto Analysis  
 Autoanalyze  
Analysis Macro

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	23.39 °C
Heat Flow	0.152 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	22.54 °C
Heater Temperature	22.02 °C
Sample Purge Flow	49.99 mL/min

#	Running Segment Description
1	↑ Equilibrate at 0.00 °C
2	↘ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↘ Ramp 10.000 °C/min to 0.00 °C
5	↑ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↘ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3


Heat Flow (mW)

Temperature (°C)

01 97.90 min

Append Apply Cancel Help

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 23.61 °C

**Experiment**  
Standard Sequence

Sequence No. 12  
Run 1:

Notes  
Operator: RICO  
Extended Text

Mass Flow Control  
Sample

Auto Analysis  
 Autoanalyze  
Analysis Macro

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	23.61 °C
Heat Flow	0.149 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	22.75 °C
Heater Temperature	22.23 °C
Sample Purge Flow	50.01 ml/min

**Autosampler Calibration - Reset auto**

Calibration | Cycle Test | Motor Test

**Reset auto**

Current Values	Down / Up	Saved Values	Home
13600	Down / Up	N/A	<input type="checkbox"/>
7085	Left / Right	N/A	<input type="checkbox"/>
0	Extend / Retract	N/A	<input type="checkbox"/>
0	Grip/Release	N/A	<input type="checkbox"/>

Tray  Arm  Optical  Electrical

Start Stop Save Close Help

01 97.90 min. Append Apply Cancel Help

Heat Flow (mW)

Temperature (°C)

Completed Standard Seg 0 in Run 1 15:47:49

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 23.71 °C

Experiment  
Standard Sequence

Sequence No. 12  
Run 1:

Notes  
Operator: RICO  
Extended Text

Mass Flow Control Settings  
Sample: #1 - Nitrogen

Auto Analysis  
 Autoanalyze  
Analysis Macro

Signal Value

Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	23.71 °C
Heat Flow	0.137 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	22.85 °C
Heater Temperature	22.34 °C
	0.00 mL/min

Calibration Set: 0) Not set

Cell Constant (Tuesday, June 07, 2011, 16:57:12)

Cell Constant: 1.0320      Onset Slope: -114.105 mW/°C

Temperature Table

	Observed Temperature °C	Correct Temperature °C
1	155.47	156.60
2	417.75	419.53
3		
4		
5		

MDSC Cp Constant (Total): 1.071  
MDSC Cp Constant (Reversing): 1.100  
Standard Cp Constant (Direct): 1.000

Reset this Set    Reset All    OK    Close    Help

Heat Flow vs. Temperature (°C) graph

01 97.90 min    Append    Apply    Cancel    Help

Completed Standard Seg 0 in Run 1 15:48:34

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The screenshot displays the QSeries software interface for a DSC instrument. The main window shows the 'Experiment' tab with 'Run 1: Complete' and a temperature of 24.67°C. A 'Calibration Report' dialog box is open, displaying the following information:

**Calibration Report**

**Last Run Calibration**

Kcell 1.0320  
Cpconst Standard: 1.000  
Calib -114.1049  
TempCal 2 pts 155.47 156.60, 417.75 419.53  
TzeroDt -90.0860 -0.0006 -1.1227  
TzeroDtz -90.0860 0.0097 3.9962  
InstCalFile Tzero: \\Ta-lab\ta\Data\DSC  
CALIBRATION\Z702\_RCS (90L\_08\_05\_2011 17\_  
34\_54\_TZR  
InstCalFile Baseline: \\Ta-lab\ta\Data\DSC  
\20110805.001  
InstCalFile Sapphire: \\Ta-lab\ta\Data\DSC  
\20110805.002  
InstCalDate Tzero 2011-08-05 Time 17:34:54  
TempRange -90.08 to 396.58 °C at 19.99 °C/min Heat  
Only  
AutoZero Delta T Offset 0.000 uV  
AutoZero Delta T0 Offset 0.000 uV  
MultiptCal 0  
Description: 0) Not set  
AutoCellConst Calibration Date 2011-06-07 16:57:12

**Multi Point Calibration Information**  
Multi Cal Index 0

Buttons: Print, Close, Help

The background interface includes a 'Signal' table with the following data:

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	24.67 °C
Heat Flow	0.137 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Sample Temperature	23.82 °C
Reference Temperature	23.30 °C
Flow	49.99 mL/min

At the bottom right, a graph shows Heat Flow (mW) vs Temperature (°C) with a scale from -10.0 to 0.0 mW and 0.90 to 2.00 °C.

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The screenshot displays the QSeries software interface for a DSC instrument. The main window is titled "QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]". The interface is divided into several sections:

- Control Panel:** Includes buttons for "Control", "Experimental", "Calibrate", and "Tools". A "Run 1: Com" indicator is visible.
- Experiment Section:** Shows "Standard Sequence" and "Sequence No. 12". A "Run 1:" indicator is present.
- Menu:** A "Tools" menu is open, listing options such as "RT Plot", "User Preferences...", "Instrument Preferences...", "Data Transfer...", "File Utility...", "Instrument Setup...", "Register as the Master Controller", "Unregister Master Controller", "Print Setup...", "Controller License...", "Instrument License...", and "Instrument Diagnostics".
- Parameters:** A "Sample" dropdown is set to "#1 - Nitrogen" and "Flow Rate" is set to "50 mL/min".
- Auto Analysis:** Includes an "Autoanalyze" checkbox and an "Analysis Macro" dropdown.
- Running Segment Description Table:**

#	Running Segment Description
1	↕ Equilibrate at 0.00 °C
2	↗ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↘ Ramp 10.000 °C/min to 0.00 °C
5	↕ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↗ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3
- Graph:** A plot of "Heat Flow (mW)" versus "Temperature (°C)". The y-axis ranges from -10.0 to 10.0, and the x-axis ranges from 0.90 to 2.00. The plot area is currently blank.
- Status Bar:** Shows "01 97.90 min" and buttons for "Append", "Apply", "Cancel", and "Help".

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 23.99°C

Help Topics...  
Getting Started...  
About Q Advantage...

### Experiment

Standard Sequence

Sequence No. 12  
Run 1:

Summary Procedure Notes

Notes

Operator: RICO

Extended Text

Mass Flow Control Settings

Sample: #1 - Nitrogen Flow Rate: 50 mL/min

Auto Analysis

Autoanalyze

Analysis Macro

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	23.99 °C
Heat Flow	0.145 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	23.14 °C
Heater Temperature	22.62 °C
Sample Purge Flow	50.02 mL/min

#	Running Segment Description
1	↕ Equilibrate at 0.00 °C
2	↗ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↘ Ramp 10.000 °C/min to 0.00 °C
5	↕ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↗ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3

01 97.90 min

Append Apply Cancel Help

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DSC Online Help - Windows Internet Explorer

Contents Index Search - Search - GO Powered By RoboHelp

TA

## Welcome to DSC Online Help

Congratulations! You are using the latest software available from TA Instruments – Waters LLC for control of our Q Series™ analyzers

Introduction

Monitoring Data

Analyzing Data

Maintaining the Instrument

Troubleshooting

Preparing the Instrument

Calibrating the Instrument

Setting Up the Experiments

See also: [Trademarks and Patents](#)

To assist you while using the online help:

As you use this help document, you can call up the menu above by clicking on the  icon in the upper right corner of each topic. You can also click the appropriate button to use the **Contents**, **Index**, **Search**, and **Print** capabilities in the help window.

### TA Instruments Web Site/MSDS Sheets

For information on the latest products, training, applications literature, software downloads, and other valuable services, visit us on the Web at <http://www.tainstruments.com>

To locate related thermal analysis and rheology Material Safety Data Sheets (MSDS) on the TA Instruments web site, click on **Support** then **Services**, and choose **Material Safety Data Sheets**.

©2001- 2006 TA Instruments—Waters LLC, 109 Lukens Drive, New Castle, DE 19720.

### Notices

The material contained in the online help for the software used to support this instrument, is believed adequate for the intended use of the instrument. If the instrument or procedures are used for purposes other than those specified herein, confirmation of their suitability must be obtained from TA Instruments. Otherwise, TA Instruments does not guarantee any results and assumes no obligation or liability. TA Instruments also reserves the right to revise this document and to make changes without notice.

TA Instruments may have patents, patent applications, trademarks, copyrights, or other intellectual property covering subject matter in this document. Except as expressly provided in written license agreement from TA Instrument, the furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property.

TA Instruments Operating Software, as well as Module, Data Analysis, and Utility Software and their associated manuals and online help, are proprietary and copyrighted by TA Instruments. Purchasers are granted a license to use these software programs on the module and controller with which they were purchased. These programs may not be duplicated by the purchaser without the prior written consent of TA Instruments. Each licensed program shall remain the exclusive property of TA Instruments, and no rights or

開始 TA Instrument Ex... QSeries - [Q400-0... Universal Analysis QSeries - [Q2000-... E - 小畫家 DSC Online Help - ... CH 3:50 PM

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

QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

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Sequence No. 12  
Run 1:

Summary Procedure Notes

Sample Information  
Sample Name: ITRI PET  
Pan Type: Tzero Aluminum  
Sample Size: 6.280 mg Pan No. 1 Ref. 1  
 Pan Mass: 0.000 mg (Sample) 0.000 mg (Reference)  
Comments

Data File Name: WTa-lab\ta\Data\DSC\LA WRENCE\TRIN\ITRI PET.001  
 Network Drive

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	22.37 °C
Heat Flow	0.173 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	21.25 °C
Heater Temperature	20.99 °C
Sample Purge Flow	50.01 mL/min

#	Running Segment Description
1	↕ Equilibrate at 0.00 °C
2	↗ Ramp 10.000 °C/min to 290.00 °C
3	✓ Mark end of cycle 1
4	↘ Ramp 10.000 °C/min to 0.00 °C
5	↕ Equilibrate at 0.00 °C
6	✓ Mark end of cycle 2
7	↗ Ramp 10.000 °C/min to 295.00 °C
8	✓ Mark end of cycle 3

01 97.90 min Append Apply Cancel Help

Heat Flow (mW)

Temperature (°C)

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 27.36 °C

### Experiment Standard Sequence

Sequence No. 12  
Run 1:

### Notes

Operator: RICO  
Extended Text:

### Mass Flow Control Settings

Sample: #1 - Nitrogen Flow Rate: 50 mL/min

### Auto Analysis

Autoanalyze  
Analysis Macro:

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	27.36 °C
Heat Flow	0.105 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	26.48 °C
Heater Temperature	25.97 °C
Sample Purge Flow	50.00 mL/min

#	Running Segment Description
1	↑ Equilibrate at 0.00 °C


### Instrument Settings [Q2000-1584 - DSC Q2000@Mfg-dsc]

DSC Name: Q2000-1584  
Serial Number: 2000 - 1584  
 Remanufactured  
Location: Mfg-dsc  
IP Address: 192 . 168 . 171 . 113  
Net Mask: 255 . 255 . 255 . 0  
Default Gateway: 0 . 0 . 0 . 0  
A/S Serial No.: 2475

Apply New Settings

< Back Finish Cancel Help


Heat Flow (mW) vs Temperature (°C)

01 97.90 min. Append Apply Cancel Help

Ready Completed Standard Seg 0 in Run 1 16:09:41

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The screenshot displays the QSeries software interface for a DSC experiment. The main window shows the 'Experiment' tab with 'Standard Sequence' selected. The 'Run 1: Complete' status is shown at the top, along with the current temperature of 26.22°C. The 'Instrument Preferences' dialog box is open, showing settings for 'MFC Purge', 'LCD Signals', 'DSC', 'Cooler', and 'Auto Sampler'. The 'Gas #1' and 'Gas #2' are both set to 'Nitrogen'. A red arrow points to the 'Gas #1' dropdown menu. The 'Stop experiment when flow rate deviates from the set value' checkbox is checked. The 'Mass Flow Control Settings' section shows 'Sample #1 - Nitrogen' and 'Flow Rate 50 mL/min'. The 'Auto Analysis' section has 'Autoanalyze' unchecked. The 'Analysis Macro' field is empty. The 'Signal' table on the right lists various parameters and their values. The 'Running Segment Description' table shows the current segment is 'Equilibrate at 0.00 °C'. The 'Heat Flow (mW)' vs 'Temperature (°C)' plot is visible at the bottom right, showing a flat line at 0.00 mW. The 'Ready' status is shown at the bottom left. The Windows taskbar at the bottom shows the system tray with the time 4:02 PM and the date 16:02:49.

QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 26.22°C

Experiment  
Standard Sequence

Sequence No. 12  
Run 1:

Summary Procedure Notes

Notes  
Operator: RICO  
Extended Text

Mass Flow Control Settings  
Sample #1 - Nitrogen Flow Rate 50 mL/min

Auto Analysis  
 Autoanalyze  
Analysis Macro

Instrument Preferences

MFC Purge LCD Signals DSC Cooler Auto Sampler

Gas #1 Nitrogen  
Gas #2 Nitrogen

Stop experiment when flow rate deviates from the set value

確定 Cancel Apply Help

Signal

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	26.22 °C
Heat Flow	0.117 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	25.37 °C
Heater Temperature	24.84 °C
Sample Purge Flow	49.99 mL/min

Running Segment Description

#	Running Segment Description
1	Equilibrate at 0.00 °C
2	min to 290.00 °C
3	min to 0.00 °C
4	min to 295.00 °C

Heat Flow (mW)

Temperature (°C)

01 97.90 min Append Apply Cancel Help

Ready

Completed Standard Seg 0 in Run 1 16:02:49

開始 TA Instrument Ex... QSeries - [Q400-0... Universal Analysis QSeries - [Q2000-... J - 小畫家 DP CH E 4:02 PM

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 26.28 °C

### Experiment

Standard Sequence

Sequence No. 12  
Run 1:

Notes

Operator: RICO

Extended Text

Mass Flow Control Settings

Sample: #1 - Nitrogen Flow Rate: 50 mL/min

Auto Analysis

Autoanalyze

Analysis Macro

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	26.28 °C
Heat Flow	0.116 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	25.42 °C
Heater Temperature	24.90 °C
Sample Purge Flow	49.99 mL/min

#	Running Segment Description
1	Equilibrate at 0.00 °C
2	Heat to 290.00 °C
3	Hold at 290.00 °C
4	Heat to 0.00 °C
5	Hold at 0.00 °C
6	Heat to 295.00 °C
7	Hold at 295.00 °C

### Instrument Preferences

MFC Purge | LCD Signals | **DSC** | Cooler | Auto Sampler

Transition Direction

Exotherm Down  **Exotherm Up**

Heat Flow Selection: **Heat Flow T4P (mW)**

Standby Temperature: 40.00 °C

Data Mode: Normal Gain

PC & Operation: Single

Enable Diagnostic Signals  AutoLid II installed

Enable Multiple Calibration Sets

確定 Cancel Apply Help

01 97.90 min Append Apply Cancel Help

Heat Flow (mW)

Temperature (°C)

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The screenshot displays the QSeries software interface for a DSC instrument. The main window shows the 'Experiment' tab with 'Standard Sequence' selected. The 'Run 1: Complete' status is shown at the top, with a temperature of 26.34 °C. The 'Instrument Preferences' dialog box is open, showing the 'Cooler' tab. The 'Cooler Selection' dropdown is set to 'RCS (90)'. The 'Activate secondary purge when lid is opened' checkbox is checked. The 'Between Runs' section has 'Leave RCS On' and 'Leave LNCS On' checked, and 'Autofill LNCS if below' set to 40%.

**Signal**

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	26.34 °C
Heat Flow	0.115 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	25.48 °C
Heater Temperature	24.96 °C
Sample Purge Flow	50.01 mL/min

**Running Segment Description**

#	Running Segment Description
1	Equilibrate at 0.00 °C
2	Heat to 290.00 °C
3	Hold at 290.00 °C
4	Cool to 0.00 °C
5	Hold at 0.00 °C
6	Heat to 295.00 °C
7	Hold at 295.00 °C

**Heat Flow (mW) vs Temperature (°C)**

Temperature (°C)	Heat Flow (mW)
0.90	0.00
1.00	0.00
1.10	0.00
1.20	0.00
1.30	0.00
1.40	0.00
1.50	0.00
1.60	0.00
1.70	0.00
1.80	0.00
1.90	0.00
2.00	0.00

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 26.45 °C

### Experiment

Standard Sequence

Sequence No. 12  
Run 1:

Summary Procedure Notes

Notes

Operator: RICO

Extended Text

Mass Flow Control Settings

Sample: #1 - Nitrogen Flow Rate: 50 mL/min

Auto Analysis

Autoanalyze

Analysis Macro

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	26.45 °C
Heat Flow	0.114 mW
Heat Capacity	0.000 mJ/°C
Set Point Temp	0.00 °C
Heater Power	0.000 W
Flange Temperature	25.60 °C
Heater Temperature	25.07 °C
Sample Purge Flow	50.01 mL/min

#	Running Segment Description
1	↑ Equilibrate at 0.00 °C
2	↘ Ramp 10.000 °C/min to 290.00 °C
3	↘ Ramp 10.000 °C/min to 0.00 °C
4	↘ Ramp 10.000 °C/min to 295.00 °C
5	↘ Ramp 10.000 °C/min to 0.00 °C

### Edit User Preferences

General | Data | Procedure | Sequence | Status

Background: TA Instruments, Inc.

Display Confirmation Message On

- Stop run
- Reject run
- Delete sequence run(s)
- Prompt to save run sequence
- Use WinUA for plot view

Buttons: 確定, Cancel, Apply, Help

01 97.90 min. Append Apply Cancel Help

Heat Flow (mW) vs Temperature (°C)

Y-axis: -10.0 to 2.0 mW  
X-axis: 0.90 to 2.00 °C

Completed Standard Seg 0 in Run 1 16:04:10

Ready

Windows Taskbar: 開始, TA Instrument Ex..., QSeries - [Q400-0..., Universal Analysis, QSeries - [Q2000-..., M - 小畫家, DP, CH E, 4:04 PM

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]

Control Experimental Calibrate Tools View Window Help

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Contents Index Search

## Message 67

### Flange over temperature. Run terminated.

**Problem:**

The flange temperature is outside the instrument's operating limit. This can occur under the following circumstances:

- If a sample starts to burn during the experiment and raises the temperature beyond the limit.
- If the cooler is not on, or hasn't been on long enough to equilibrate (e.g., flange temperature is too high).
- If the cooler is not operating correctly (e.g., has run out of liquid nitrogen).

**Solution:**

1. Try running the sample in an inert atmosphere and/or to a lower temperature.
2. Turn the cooler on at least 20 minutes prior to starting the run.
3. Check the operation of the cooler.

Call TA Instruments for service, if the problem persists.

General Instrument Messages

- ? What are General Instrument Mess
- ? Message 10
- ? Message 21
- ? Message 60
- ? Message 62
- ? Message 63
- ? Message 66
- ? Message 67
- ? Message 73
- ? Message 75
- ? Message 78
- ? Message 80
- ? Message 81
- ? Message 91
- ? Message 96
- ? Message 97
- ? Message 98
- ? Message 99
- ? Message 601
- ? Message 603
- ? Message 605
- ? Message 606
- ? Message 674
- ? Message 675
- ? Message 676
- ? Message 677
- ? Message 678
- ? Message 679
- ? Message 680
- ? Message 681
- ? Message 682
- ? Message 684

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QSeries - [Q2000-1584 - DSC Q2000@Mfg-dsc]


Control Experimental Calibrate Tools View Window Help

Run 1: Complete Temp: 26.82°C

### Experiment

#### Standard Sequence

Sequence No. 12  
Run 1:

This wizard helps you set up and start typical DSC experiments. The types of experiments you can run depend on the Q Series DSC and cooling accessory you are using. Options which do not apply for your DSC system (DSC and cooling accessory) are automatically "grayed out" so they cannot be inadvertently selected.

Conventional DSC  
These experiments measure heat flow into or out of a material as a function of temperature or time.

Modulated DSC ® (MDSC ®)  
These experiments provide the same heat flow information as conventional DSC, but also separate that heat flow into heat capacity related (reversing) and kinetic (nonreversing) components.

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Purge Gas

✓ Be sure that your purge gas is connected and properly regulated. Purge Gas is recommended for all DSC experiments.

Cooling Accessory

✓ The Q Series DSC cells generally require a cooling accessory to be connected regardless of the type of experiment being run. The QCA is the only exception. It is removed before initiating heating ramp experiments.

Loading the Sample

✓ Position the sample pan into the cell. [Typically, an empty cell (e.g., no pans) is used for this calibration.] Position a reference pan, of the same type, into the cell (back position). Cover the cell with the appropriate lids.

Append Run Start Run Finish Cancel Help

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The screenshot displays the QSeries software interface for a DSC Q10. The main window shows a graph of Heat Flow (mW) versus temperature, with a 'Run 1: Standby Temp: 40.85°C' label. An 'Open Calibration File' dialog box is open, showing a file list with the following items:


- Baseline 980219
- Dry Cell condition 980219
- Replication

The 'Open File' button in the dialog is highlighted with a red box. A red text annotation '選取欲分析的檔案' (Select the file to be analyzed) is placed over the file list. The software interface includes a menu bar (Control, Experimental, Analysis, Tools, View, Window, Help) and a toolbar with various icons. The status bar at the bottom shows 'Ready', 'Stand by', 'Calibration', 'Seg 0 in Run 1', and the time '14:00:56'.

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
## Conventional DSC Experiments [Q10-0300 - DSC Q10@Mfg-dsc]

Select the type of DSC calibration experiment you want to perform:

- Baseline  
This calibration compensates for subtle differences between the reference and sample thermocouples. The baseline calibration is based on heating an empty cell through the same temperature range that will be used in subsequent experiments. Baseline calibration is recommended any time the heating/cooling rate, purge gas, or cooling accessory is changed.
- Cell Constant / Temperature  
Cell Constant is a calibration factor used to adjust for subtle differences in the calorimetric response of a DSC cell. Temperature calibration ensures that the sample thermocouple reading is correct under the experimental conditions chosen. Both of these calibrations are performed based on the melting peak of a standard metal such as indium. These calibrations should be performed any time that the heating/cooling rate, purge gas, cooling accessory, or pan type is changed.

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
Experimental Parameters: Cell Constant [Q10-0300 - DSC Q10@Mfg-dsc]

Cell constant and temperature calibrations involve heating a high purity, metal standard (e.g., indium) through its melting peak using the same ramp rate, purge gas and pan type that will be used for subsequent experiments.

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< Back Next > Cancel Help

Signal

- Method Time
- Segment Time
- Remaining Run Time
- Temperature
- Delta T
- Heat Flow
- Set Point Temp
- Heater Power

Enter Desired Data File Name

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- Baseline 980219
- Dry Cell condition 980219
- Replication

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
**Purge Gas**

✓ Be sure that your purge gas is connected and properly regulated. Purge Gas is recommended for all DSC experiments.

**Cooling Accessory**

✓ The Q Series DSC cells generally require a cooling accessory to be connected regardless of the type of experiment being run. The QCA is the only exception. It is removed before initiating heating ramp experiments.

**Loading the Sample**


✓ Position the sample pan into the cell. [Typically, an empty cell (e.g., no pans) is used for this calibration.] Position a reference pan, of the same type, into the cell (back position). Cover the cell with the appropriate lids.

Append Run **Start Run** Finish Cancel Help

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The screenshot displays the QSeries software interface for DSC Q10 calibration. The main window shows a graph with 'Heat Flow (mW)' on the x-axis (0.90 to 2.00) and temperature on the y-axis (0.90 to 2.00). A status bar at the top indicates 'Run 1: Complete Temp: 40.82°C'. An 'Open Calibration File' dialog box is open, showing a search path of '990219' and a list of files including 'Baseline 990219', 'Dry Cell condition 990219', 'Indrum 990219', 'Indrum 990219', and 'Replication'. The 'Indrum 990219' file is selected. The dialog box also shows the file name 'Indrum 990219' and the file type 'Calibration Files (\*.\*)'. The 'Open File' button in the bottom-left corner of the main window is circled in red.

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DSC calibration should be performed using the same conditions (purge gas, sample pan type, and cooling accessory) that will be used in subsequent experiments on your sample materials. In addition, the calibration sequence performed depends on the type of heat flow that will be stored in your subsequent experiments. Selection of the highest level heat flow for your Q Series DSC is recommended.

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DSC calibration should be performed using the same conditions (purge gas, sample pan type, and cooling accessory) that will be used in subsequent experiments on your sample materials. In addition, the calibration sequence performed depends on the type of heat flow that will be used in your experiments. Selection of the highest level heat flow for your Q Series DSC is recommended.

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This wizard helps you setup and start experiments used to calibrate for cell resistance & capacitance differences, heat flow (cell) constant, and temperature. Select from one of the calibration options.

If your instrument is equipped with an autosampler, it will be disabled during this calibration.

Cell Resistance & Capacitance, Cell Constant, and Temperature Calibration:

Cell resistance and capacitance calibration compensates for subtle differences in thermal resistance and capacitance between the reference and sample platforms in the DSC sensor. The calibration is based on two experiments - one run with an empty cell, and the second run with equal weight sapphire disks on the sample and reference platforms. This calibration is recommended any time the purge gas or cooling accessory used in subsequent experiments are changed.

Cell Constant & Temperature Calibration

Cell constant is a calibration factor used to adjust for subtle differences in the calorimetric response of a DSC cell. Temperature calibration ensures that the sample thermocouple reading is correct under the experimental conditions chosen. Both of these calibrations are performed on the melting peak of a standard metal such as indium. These calibrations should be performed any time that the heating/cooling rate, purge gas, or cooling accessory are changed.

Load Saved Cell Resistance & Capacitance Calibration File to Instrument

This option allows you to send a saved cell resistance and capacitance (Tzero) calibration results file to the instrument.

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Cell resistance & capacitance calibration involves two experiments. The first experiment involves heating an empty cell through the widest temperature range that will be used in subsequent experiments.

Conduct experiment

1. Remove any pans from cell and cover cell.
2. Enter the desired test parameters:

Calibration range: Tzero Heat only cal

Lower temperature: -90.00 °C

Upper temperature: 400.00 °C

Ramp rate: 20.000 °C/min

Operator: Rico

Comment: Tzero Calibration - Baseline Run

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
The second calibration experiment involves heating the cell with two equal weight sapphire disks on the sample and reference platforms through the entire temperature range under the same conditions as the previous experiment. The sapphire disks are placed directly on the cell platforms (i.e., no sample pans are used).

Conduct experiment

1. Place sapphire 1 on sample platform in cell.

2. Enter weight of sapphire 1:  mg **1.輸入透明寶石重量**

3. Place sapphire 2 on reference platform in cell.

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Run 1: Complete Temp: 39.93°C

Cell constant and temperature calibrations involve heating a high purity, metal standard (e.g., indium) through its melting peak using the same ramp rate, purge gas and pan type that will be used for subsequent experiments. For the most accurate results, calibration should be done with the highest level heat flow signal available in the instrument.

1. Remove sapphire samples from cell.
2. Select calibration standard for cell constant calibration

Standard: Indium Weight: 95.160 mg

Pan Type: Tin

Lead

Zinc

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
3. Place sample pan containing standard on front platform in cell.
4. Place reference pan on back platform in cell.

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Cell constant and temperature calibrations involve heating a high purity, metal standard (e.g., indium) through its melting peak using the same ramp rate, purge gas and pan type that will be used for subsequent experiments. For the most accurate results, calibration should be done with the highest level heat flow signal available in the instrument.

1. Remove sapphire samples from cell.
2. Select calibration standard for cell constant calibration

Standard: Indium Weight: 1.45 mg

Pan Type: Tzero Aluminum

3. Place sample pan on front platform in cell.
4. Place reference pan on back platform in cell.

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Conduct experiment

1. Enter the desired test parameters for cell constant and temperature calibration.

Start Temperature:  Use current  Premelt  
100.00 °C 1. 設定起始溫度

Heating Rate: 10.000 °C/min 2. 設定升溫速率

Final Temperature: 180.00 °C 3. 設定結束溫度

Data File Name: \\Lotus\sta\Data\DSC\980217.003 4. 設定存檔路徑

Network Drive

Operator: Rico

Purge Gas: #1 - Nitrogen Flow Rate: 50 mL/min 5. 設定氮氣流量

2. Cover Cell.

Enter existing Cell Constant data file:

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

Analysis complete. Calibration results successfully saved to instrument.

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