# Discovery TGA 55xx Series 操作訓練

3

**2018.05.31**

**Allen Chang**

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# Discovery TGA 55xx Series 操作訓練說明

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

**WARNING: Oxygen can be used as a purge gas in the TGA. However, if you use oxygen as a purge gas, you must make sure the furnace is cleaned of hydrocarbons that could combust.**

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




The screenshot displays the TGA control software interface. The main window is titled "Experiments" and contains several sections:

- File Manager:** Located on the left, it shows a list of experiments under "Experiments", "Incomplete Queue (0)", "Running Queue (1)", and "Design View (1)".
- Experiments View:** The central area shows the details of a selected experiment, including "Notes", "File Name", "Procedure", "Segments", and "Advanced" settings.
- Control Panel:** Located on the right, it displays the current status of the instrument, including "Idle 26.85 °C", a table of signals, and various control parameters.

A red box highlights the "Experiments" view window, and a red arrow points to it with the text "Experiments視窗".

Signal	Value	Units
Method Time	0.00	min
Remaining Segment Time	0.00	min
Remaining Time	0.00	min
Temperature	26.85	°C
Weight	843.581	mg
Temperature Difference	0.176369	°C
Sample Purge	60.00	mL/min
Balance Purge	38.31	mL/min
Set Point Temperature	26.64	°C
Power Delivered	0.0001	W
Activation Energy (kJ/mol)		
Power Requested	0.0000	W
Ln(Rate Ratio)		
Log Pre-Exponential Factor		Time

Class	Value	Units
Gas 1 Nitrogen	60.00	mL/min
Balance Flow	38.31	mL/min
Set Sample Flow	25.00	mL/min
Gas 1 Nitrogen		

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The screenshot displays the TGA control software interface. The main window is titled "Experiments" and shows a procedure editor. The procedure is named "Run 1 - (Plan 1 - Custom)" and consists of two segments:

No.	Description
1	Equilibrate 50.00 °C
2	Ramp 10 °C/min to 500.00 °C

The right-hand side of the interface features a "Control Panel" showing the current temperature as 26.85 °C. Below this, there is a table of signals and their values:

Signal	Value	Units
Method Time	0.00 min	
Remaining Segment Time	0.00 min	
Remaining Time	0.00 min	
Temperature	26.85 °C	
Weight	8.43581 mg	
Temperature Difference	0.176369 °C	
Sample Purge	60.00 mL/min	
Balance Purge	38.31 mL/min	
Set Point Temperature	26.64 °C	
Power Delivered	0.0001 W	
Activation Energy (kJ/mol)		
Power Requested	0.0000 W	
Ln(Rate Ratio)		
Log Pre-Exponential Factor		1 min

A red box highlights the "Experiments" menu item in the bottom-left corner, with a red arrow pointing to it from a red box containing the text "功能模式切換視窗".

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The screenshot displays the TGA control software interface. The main window is titled "Experiments" and shows a "Design View" for a test run. The "Segments" table is visible, listing two segments: "1 Equilibrate 50.00 °C" and "2 Ramp 10 °C/min to 500.00 °C". The "Advanced" section includes checkboxes for "Start Experiment After Weight Stabilization" and "Enable Air Cool", with "Air Cool For" set to 5.0 min. A red box highlights the "Experiments" window and the "Control Panel" on the right. The "Control Panel" shows the current temperature as "Idle 26.85 °C" and a table of signals.

Signal	Value	Units
Method Time	0.00 min	
Remaining Segment Time	0.00 min	
Remaining Time	0.00 min	
Temperature	26.85 °C	
Weight	8.43581 mg	
Temperature Difference	0.176369 °C	
Sample Purge	60.00 mL/min	
Balance Purge	38.31 mL/min	
Set Point Temperature	26.64 °C	
Power Delivered	0.0001 W	
Activation Energy (kJ/mol)		
Power Requested	0.0000 W	
Ln(Rate Ratio)		
Log Pre-Exponential Factor		1 min

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The screenshot shows the TA Instruments TRIOS v4.1.0.31.739 software interface. The main window displays the 'Experiments' section with a 'Design View' showing a procedure with two segments: '1 Equilibrate 50.00 °C' and '2 Ramp 10 °C/min to 500.00 °C'. A 'Control Panel' window is open on the right, showing the current status 'Idle 26.85 °C' and a table of signals.

Signals [14 of 14]	Value	Units
Method Time	0.00 min	
Remaining Segment Time	0.00 min	
Remaining Time	0.00 min	
Temperature	26.85 °C	
Weight	8.43 581 mg	
Temperature Difference	0.176369 °C	
Sample Purge	60.00 mL/min	
Balance Purge	38.31 mL/min	
Set Point Temperature	26.64 °C	
Power Delivered	0.0001 W	
Activation Energy (kJ/mol)		
Power Requested	0.0000 W	
Ln(Rate Ratio)		
Log Pre-Exponential Factor		1 min

A red box highlights the 'Control Panel' window, and a red arrow points from a central label to it.

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The screenshot displays the TA Instruments TGA control software interface. The main window is titled "Experiments" and shows a procedure configuration for a "Custom" test. The procedure consists of two segments: 1. Equilibrate 50.00 °C and 2. Ramp 10 °C/min to 500.00 °C. The "Advanced" section includes options for "Start Experiment After Weight Stabilization" and "Enable Air Cool" with a 5.0 min air cool factor.

On the right side, the "Control Panel" shows the current temperature at 26.85 °C and a list of signals. A red box highlights the "General" settings section, which includes:


- Gas 1 Nitrogen: 60.00 mL/min
- Balance Flow: 38.31 mL/min
- Set Sample Flow: 25.00 mL/min
- Gas 1: Nitrogen

A red arrow points from the Chinese text "儀器一般功能設定視窗" (Instrument General Function Setting Window) to the highlighted "General" settings section.

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





The screenshot displays the software interface for TGA 5500, specifically the Calibration Data section. The interface includes a navigation bar with tabs for Calibration Data, Temperature Calibration, DTA Signal Setup, and Weight Calibration. The main content area is divided into three sections: Temperature Calibration (Not Applied), DTA Signal (Not Applied), and Weight Calibration. The Weight Calibration section contains a table with Autocalibration Pan Information and a table of calibration data.

**Calibration Data**

Temperature Calibration (Not Applied)

DTA Signal (Not Applied)

**Weight Calibration**

Autocalibration Pan Information

Calibration Fixture 1 Pan Number	11
Calibration Fixture 2 Pan Number	12
Calibration Fixture 3 Pan Number	13

Weight	Offset	Gain
443.436 mg	42.2906 mg	0.0000 mg/count
1243.055 mg	639.4948 mg	0.0001 mg/count

Information (4013) 4/13/2017 10:29:20 AM

Weight calibration has successfully completed.

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The screenshot displays the software interface for TGA 55X. The main window is titled "Calibration" and contains several tabs: "Calibration Data", "Temperature Calibration", "DTA Signal Setup", and "Weight Calibration". The "Calibration Data" tab is active, showing sections for "Temperature Calibration (Not Applied)", "DTA Signal (Not Applied)", and "Weight Calibration".

An information dialog box is overlaid on the screen, titled "Information (4013) 4/13/2017 10:29:20 AM". It contains an information icon and the text "Weight calibration has successfully completed." with a "Close" button.

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





The screenshot displays the 'Temperature Calibration' software interface. The main window shows the following details:

- Heating Rate:** 10 °C/min
- Calibration Values Table:**

Standard	Material Type	Reference	Measured	Offset
None	Cure Point	200.0	NaN	NaN

- Calibration Input Data Table:**

Standard	Material Type	Reference	Measured	Apply Offset	Value
None	Cure Point	200.0	NaN	Off	Go to

The right-hand 'Control Panel' displays the current status as 'Idle 222.29 °C' and lists various signals and their units, including 'Flowing Segment Time', 'Flowing Time', 'Temperature', 'Weight', 'Temperature Difference', 'Sample Purge', 'Balance Purge', 'Set Point Temperature', 'Power Delivered', 'Autation Every 60 Min', 'Power Resistor', 'LoFlow Ratio', and 'Loop Pre Exponential Factor'. The bottom 'General' section shows 'Gas 1 Nitrogen' and 'Release Flow' set to 80.00 mL/min.

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
The screenshot displays the 'Control Panel' window with the following data:

Signals (14 of 14)	Value	Units
Method Time	0.00	min
Remaining Segment Time	0.00	min
Remaining Time	0.00	min
Temperature	26.85	°C
Weight	#-43.581	mg
Temperature Difference	0.175368	°C
Sample Purge	60.00	mL/min
Balance Purge	38.31	mL/min
Set Point Temperature	26.64	°C
Power Delivered	0.0001	W
Activation Energy (kJ/mol)		
Power Requested	0.0000	W
Ln(Rate Ratio)		
Log Pre-Exponential Factor		1/min

Below the table, the 'General' section shows the following settings:

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The bottom of the interface shows a status bar with 'General' and 'Auto Sampler' tabs, and a system tray with the time '2:03 PM' and date '4/12/2017'.

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